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PLEASE READ CAREFULLY

# AntiCruZ

(Indirect immunofluorescence assay for the detection of *Trypanosoma cruzi* antibodies)

Only for professional use

### Application

In this indirect immunofluorescence assay slides with coated *Trypanosoma cruzi* infected rat heart are used to detect antibodies from human sera of patients with potential of Chagas' disease.

### Principle

In a first step the diluted serum will be incubated on a slide with coated heart muscle cells. If antibodies against *T. cruzi* antigen in the patient sample are present, they will bind to Pseudocytan of the heart muscle. In a washing step unbound serum components will be removed. By incubation with fluorescein isothiocyanate (FITC) conjugated antibodies which are directed against immunoglobulin of class G, A, and M, the bound antibodies of the serum sample can be specifically labeled. After a further washing step and embedding the tissue section it can be examined under a fluorescence microscope.

### Storage and stability

 light protected 2-8°C (fridge)

 10 months from date of manufacture

### Content

- Slides with coated antigen fields
- Serum positive control, ready to use
- Serum negativ control, ready to use
- FITCIg(G, A, M)-conjugate with Evans-Blue, ready to use
- Mounting medium, ready to use
- Blotting paper
- Dilution and washing buffer

### Required utensils and materials

- Distilled water
- Gloves
- Staining cuvettes and slide holder
- Pencil
- Fluorescence microscope
- Test tubes for serum dilution
- Humidity chamber
- Coverslips 24 x 60

### Preparation

Heat the required number of slides, the reagents and the sera to room temperature. Prepare the buffer solution by completely dissolving the entire contents of the bag with the PBS buffer salts in 1000 ml deionized water. Remove slides from the packaging. According to the protocol label them on the mat part of the slices with a pencil and put them in the humid chamber. Prepare the sera in a ration 1:10 with Diluent (for example, 5 µl serum in 45 µl dilution buffer).

### Procedure (see also scheme next page)

**Throughout the implementation the tissue slices must not dry out!**

**Preincubation:** Apply 10 µl of sample buffer to each to be examined fields on the slide, incubate 10 min. and remove the buffer with absorbent paper.

1. Apply a drop of the positive and negative control sera, 15 µl of diluted patient samples to each antigen field on the slide(s). Pay attention to the complete coverage of the tissue section.
2. Incubate the slide(s) 60 min. at 37°C in a humid chamber. Then aspirate with blotting paper the sera without touching the tissue sections.
3. Place the slide(s) into the slide holder and immerse briefly once in the staining buffer filled cuvette. Change the buffer, and then wash the slide(s) twice for 5 min. Wash and take the slide(s) out of the holder.
4. Dry the slide(s) between the antigen fields and apply 25-50 µl of FITC solution (ready to use) to each antigen fields (ensure complete coverage).
5. Incubate the slide(s) 30 min. at 37°C in a humid chamber. Then wash the slide(s) as stated in step 3.
6. Remove the slide(s) from the slide holder and remove excess buffer solution. Add 2 drops of mounting medium to each slide and place the cover glass without damaging the tissue section and without any air bubbles.
7. The evaluation is performed at 400x magnification under a fluorescence microscope.

## Goals of evaluation

1. The result is negative if no or a weak fluorescent background staining is visible. Because of counterstaining the tissue is stained weakly reddish.
2. Investigate the slides within an hour. If the slides are kept refrigerated in a moist chamber, they can be investigated within 24 hours. Let refrigerated slides reach room temperature before you perform the investigation.
3. Examine the slides at 400x total magnification.
4. Drying of slides can disturb the most peripherally situated organisms in the well, which is why you should disregard the reaction of these organisms.
5. The reaction should only be regarded as positive if the majority of the parasites (>80%) have a significant homogeneous surface fluorescence.
6. The intensity of the color reaction of individual pseudocysts may vary. Nevertheless, the degree of coloration is based on the general appearance of the antigen.

If a quantitative statement for a positive result is required, the serum should be titrated.


## Safety information

Due to the reagents contain sodium azide as a preservative, skin contact should be avoided. The sera included in the kit were tested as HBsAg and HIV antibody negative. Nevertheless, they should be handled as potentially infectious and the usual laboratory safety measures should be carried out.

## References

1. Gross, R., Schölmerich P. Lehrbuch der Inneren Medizin, 8. Auflage, Schattauer Verlag Stuttgart – New York 1994
2. Wick G., Baudner S., Herzog F. Immunofluoreszenz Die Medizinische Verlagsgesellschaft, Marburg/Lahn 1976
3. Friemel H. Immunologische Arbeitsmethoden Gustav Fischer Verlag Stuttgart 1984
4. Bittencourt, AL. Congenital Chagas disease. Am J Dis Child 1976 130:97-103.
5. Koberle F. Chagas' disease and Chagas' syndromes: The pathology of American trypanosomiasis In: Dawes B. (Ed) Advances in Parasitology 6:63-116 Academic Press, New York
6. Hoff R. et al. *Trypanosoma cruzi* in the cerebrospinal fluid during the acute stage of Chagas' disease. N Eng J Med 1978 298:604.
7. Cossio PM et al. Chagasic cardiomyopathy. Am J Pathol 1977 86:533
8. Martini-Campos JV, Tafuri WL. Chagas' enteropathy. Gut 1973 14:910. Northeast Brazil. Am J Trop Med Hyg. 1978 27:1123-1127.
9. Wkly Epidem Rec WHO 1985 60:37-44.
10. Cerisola JA. Immunodiagnosis of Chagas' disease: hemagglutination and immunofluorescence tests. J Parasitol. 1970 56:409-10

## Scheme (read also chapter „procedure“)

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