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For Professional use Only

Application

The determination of the content of fructose in the seminal fluid is a biochemical marker for the secretory function of the seminal vesicle (seminal vesicle). It is also an additional diagnostic value in terms of transport abnormalities in sperm or closures of the seminal tract (high shutter near the urethra or vas deferens aplasia).

Principle

Under the influence of concentrated hydrochloric acid (low pH), combined heat, fructose forms in the presence of indole a color complex, which can be determined spectrophotometrically.

Storage and stability

- ☒ 2-8°C (light protected)
- ☒ 24 months from date of manufacture

Content

▪ Reagent 1 (TCA)	10 ml
▪ Reagent 2 (Fructose)	10 ml
▪ Reagent 3 (Indole)	5 ml
▪ Reagent 4 (HCl)	2 x 10 ml
▪ Reagent 5 (NaOH)	5 ml
▪ Micro titer plate	

Required utensils (not included)

- Distilled water
- Gloves
- Native ejaculate or seminal plasma (100 µl)
- Slides
- Paper towels
- Pipettes and tips (10-100 and 100-1000 µl)
- Test tubes (1.5 or 2 ml)
- Test tube holder
- Centrifuge
- Micro-plate reader
- Water bath or incubator 37°C
- Vortex / mixer

Examination of fructose standard curve

		Dist. H ₂ O	Reagent 2	Fructose- concentration
1. Blank	V ₀	100 µl	0 µl	0 mg/ml
2. Dilution	V ₁	80 µl	20 µl	1 mg/ml
3. Dilution	V ₂	60 µl	40 µl	2 mg/ml
4. Dilution	V ₃	40 µl	60 µl	3 mg/ml
5. Dilution	V ₄	20 µl	80 µl	4 mg/ml
6. Dilution	V ₅	0 µl	100 µl	5 mg/ml

It is recommended to re-create the standard curve in each experiment.

Procedure

1. Spin ejaculate 10 min at 1000xg.
2. Transfer 100 µl of the supernatant (seminal plasma) in a new tube
3. Add 100 µl of reagent 1 and incubate 15 min at room temperature.
4. Centrifuge then 15 min at 1000xg.
5. Transfer carefully the supernatant in a new reagent (sample is free of protein) and store in a cooled tube holder.
6. Pipette 25 µl of this sample and of standard dilution series in new tubes.
7. Pipette 25 µl of reagent 3 to the sample and to the dilution series and mix.
8. Pipette 200 µl of reagent 4 to all test tubes.
9. Close the tubes and incubate 30 min in a water bath or 60 min in an incubator
10. Stop the reaction with 200 µl of reagent 5.
11. Centrifuge all tubes for 1 min at 1000xg.
12. Pipette 200 µl of the supernatant (see scheme below) to the microtiter plate and measure the OD values (optical density, absorbance) at 470 nm in a micro-plate reader.

	1	2	3	4	5	...	12
A	V_0	V_0	P_1	P_1			
B	V_1	V_1	P_2	P_2			
C	V_2	V_2	P_3	P_3			
D	V_3	V_3					
E	V_4	V_4					
F	V_5	V_5					
G							
H							

Standard dilution: A1 and A2 to F1 and F2

Evaluation goal: To determine the fructose amount in total ejaculate.


The average value is 13 µMol of fructose (2.4 mg) per ejaculate. If the amount is less than 13 µMol probably a blockage of the vas deferens, a congenital absence of the vas deferens, a partial reflux of ejaculate, an underactive glands within the seminal vesicle or an androgenic deficiency impairment is possible (WHO, 2010).

I. OD of the standard curve and samples:

	1	2	3	4	5	...	12
A	0.07	0.06	0.29*	0.31*			
B	0.10	0.09	0.20	0.21			
C	0.12	0.11	0.22*	0.38*			
D	0.16	0.15					
E	0.19	0.18					
F	0.22	0.22					
G							
H							

Creating of a fructose standard curve:

The fructose standard curve is created using a spreadsheet program (e.g. Microsoft Excel®). On the x-axis, the fructose concentrations (mg/ml) and on the y-axis the OD values were plotted. This also allows the calculation of the coefficient of determination (R^2), which is used to determine the working precision (1 to 99% and from 0.01 to 0.99) in the preparation of the standard curve; in the case of lower values than 0.99 the standard curve should be re-created.

*ODs of samples which are not within the range of the standard curve cannot be used for concentration determination. This must then be repeated with adjusted sample volumes.

II. Calculation of sample concentration:

The concentration of the sample can now be read or determined by using the slope function. The general formula of the standard curve is:

$$y = ax + c$$

y = OD sample

x = concentration of sample used (mg/ml)

a = slope function

c = linear constant of the curve

III. Calculation of the dilution factor:

Total volume (450 µl): (12.5 µl seminal plasma + 12.5 µl reagent 1 + 25 µl reagent 3 + 200 µl reagent 4 + 200 µl reagent 5) / 12.5 µl = 36

IV. We are looking for the total amount of fructose per ml of ejaculate. The calculation is performed in three steps:

1. Calculation of the concentration x of the sample used, derived from the general formula of the absorbance:

$$x \text{ (mg / ml)} = (y - c) / a$$

2. Calculation of the concentration of the undiluted sample

$$x \cdot \text{dilution factor}$$

3. To determine the total amount of fructose, the result is to be multiplied by the total volume (V) in ml of seminal plasma.

$$x \cdot \text{dilution factor} \cdot V$$

Note: The standard curve is linear up to 0.5 mg/ml; lower fructose concentrations cannot be accurately determined with this test.

Example

We are looking for the fructose amount of the discharged semen of patients P1 (3 ml) and P2 (5 ml). The calculation is performed in three steps:

1. Calculation of the amount of fructose in samples

$$x \text{ (mg/ml)} = (y - c) / a$$
$$P_2 \text{ (mg/ml)} = (0.205 - 0.065) / 0.03 \text{ ml/mg} =$$

4.66 mg/ml

2. Calculation of the amount of fructose in diluted samples

$$x / \text{ml} \cdot \text{dilution factor} = 4.66 \text{ mg/ml} \cdot 36 =$$

167.76 mg/ml

3. Calculation of the amount of fructose in ejaculate volume

amount of fructose (not diluted) • ejaculate volume

$$167.76 \text{ mg/ml} \cdot 5 \text{ ml} = \mathbf{838.8 \text{ mg}}$$

Results: The sample of patient 2 is noticeable.

Example 1: OD value above the highest standard

	old	new
Volume sample	25 µl	12.5 µl
Reagent 3	200 µl	237.5 µl

After measuring, the result is to be multiplied by a factor of 73. Total volume (450 µl): (6.25 µl seminal plasma + 6.25 µl reagent 1 + 25 µl reagent 3 + 212.5 µl reagent 4 + 200 µl reagent 5) / 6.125 = **73**

Example 2: OD below the lowest standard

	old	new
Volume sample	25 µl	50 µl
Reagent 3	200 µl	175 µl

After measuring, the result is to be multiplied by a factor of 18. Total volume (450 µl): (25 µl seminal plasma + 25 µl reagent 1 + 25 µl reagent 3 + 175 µl reagent 4 + 200 µl reagent 5) / 25 = **18**

Safety information / precautions (Please also consult SDS)

- All semen samples should be considered potentially infectious.
- Handle with all samples like HIV or hepatitis infected material.
- When working with samples and reagents wear always protective clothing (gloves, gowns, eye / face protection).
- Reagent 1 contains Trichloroacetic acid solution. It causes severe burns, is very toxic to aquatic organisms and may cause long-lasting adverse effects on the aquatic habitat. In case of eye contact, immediately rinse with plenty of water and seek medical advice. In

case of accident or if you feel unwell seek immediate medical attention.

- Reagent 3 (Indole) is harmful if swallowed and highly flammable. Avoid skin and eye contact.
- Reagent 4 (HCl). It causes burns on contact and irritation of the skin and respiratory tract irritation. In case of eye contact, immediately rinse with plenty of water and seek medical advice. In case of accident or if you feel unwell seek immediate medical attention.
- Reagent 5 (NaOH) causes burns. In case of eye contact, immediately rinse with plenty of water and seek medical advice. In case of accident or if you feel unwell seek immediate medical attention.

Never put water to concentrated acid (danger of explosion!).

References

1. Cooper TG et al. (1991) Variations in semen parameters from fathers. *Human Reproduction*, 6:859-866
2. Daudin M et al. (2000) Congenital bilateral absence of the vas deferens: clinical characteristics, biological parameters, cystic fibrosis transmembrane conductance regulator gene mutations, and implications for genetic counseling. *Fertility and Sterility*, 74:1164-1174
3. Eckardstein S et al. (2000) Seminal plasma characteristics as indicators of cystic fibrosis transmembrane conductance regulator (CFTR) gene mutations in men with obstructive azoospermia. *Fertility and Sterility*, 73:1226-1231
4. Karvonen MJ, Malm M (1955) Colorimetric determination of fructose with indol. *Scandinavian journal of clinical and Laboratory Investigation*, 7; 305-307
5. De la Taille A et al. (1998) Correlation of genitourinary abnormalities, spermiogram and CFTR genotype in patients with bilateral agenesis of the vas deferens. *Progress in Urology*, 8:370-376
6. WHO Press, (2010) Laboratory manual for the examination and processing of human semen. 5th. Edition

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Referenzen

1. Cooper TG et al. (1991) Variations in semen parameters from fathers. *Human Reproduction*, 6:859-866
2. Daudin M et al. (2000) Congenital bilateral absence of the vas deferens: clinical characteristics, biological parameters, cystic fibrosis transmembrane conductance regulator gene mutations, and implications for genetic counseling. *Fertility and Sterility*, 74:1164-1174
3. Eckardstein S et al. (2000) Seminal plasma characteristics as indicators of cystic fibrosis transmembrane conductance regulator (CFTR) gene mutations in men with obstructive azoospermia. *Fertility and Sterility*, 73:1226-1231
4. Karvonen MJ, Malm M (1955) Colorimetric determination of fructose with indol. *Scandinavian journal of clinical and Laboratory Investigation*, 7: 305-307
5. de la Taille A et al. (1998) Correlation of genitourinary abnormalities, spermogram and CFTR genotype in patients with bilateral agenesis of the vas deferens. *Progress in Urology*, 8:370-376
6. WHO Press, (2010) Laboratory manual for the examination and processing of human semen. 5th. Edition