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Professional use only

Application

RapidStain is a quick-staining-method to assess the morphology of smear. This method is composed of a staining kit which allows differential staining of the smear parts due to their different basophilic, eosinophilic and neutrophilic properties.

Principle

The smears are fixed in the first reagent. Here, the succession staining is used, which means, two dyes are used one after the other and it comes to differentiated staining of different tissues with individual dyes.

Storage and stability

15-25°C

36 months from date of manufacture

Content

- Reagent 1 1x 50 or 250 ml
- Reagent 2 1x 50 or 250 ml
- Reagent 3 1x 50 or 250 ml

Necessary material

- Sample (5-10 µl)
- 4 staining cuvettes, or coplin jars
- Gloves
- Tweezers
- Paper towels
- Slides
- Slides rack (if more than five slides are to dye)
- Tap or distilled water
- Immersion oil
- Microscope

Procedure (see also scheme next page)

1. Apply 5-10 µl smears per slide, grease the smears with a coverslip and let dry.

2. Fill the staining cuvettes with reagent 1, reagent 2 and reagent 3
3. Immerse the slide 5x in reagent 1
4. Flow away reagent 1 from the slide
5. Immerse the slide 5x in reagent 2
6. Flow away reagent 2 from the slide
7. Immerse the slide 5x in reagent 3
8. Flow away reagent 3 from the slide
9. Immerse the slide 20x in tap or distilled water
10. Dry the slide completely on air
11. Observe the slide under a microscope with an objective of 1000x and oil immersion

Evaluation

Reading and interpretation of results are for use by qualified professional staff. After the evaluation the immersion oil can be gently removed from the slide with a paper towel. Then the slide can be immersed in reagent 1 for 5 min, dried and stored. It is also possible to produce preparations with a coverslip and glue for long term storage.

Structure	Color
Blood smear (Haematology)	
Erythrocytes	reddish
Platelets	purple nucleus
Granulocytes	light purple granule
Granulocytes	brick-red to reddish
Eosinophils	brown granule
Monocyte	plasma light blue
Blood smear (Malaria see Fig. 1)	
Erythrocytes with one or more granules	
Sperm morphology	
Nucleus	violet
Acrosome	pink
Middle piece	reddish to pink

Warnings and precautions

- All smear samples should be considered potentially infectious.
- Handle with all samples like HIV or hepatitis infected material.
- When working with samples and reagents wear always protective clothing (gloves, lab coat, eye / face protection).
- Reagent 1 is containing methanol: toxic by inhalation, skin contact or ingestion. May cause organ damage. There is a risk of irreversible damage.
- All other ingredients are not classified as toxic



Fig. 1: Erythrocytes with *Plasmodium berghei* in blood cells (©ANTIGENES)

Scheme (see also chapter “method”)

References

- Andersen AG et al.** (2000) High frequency of sub-optimal semen quality in an unselected population of young men. *Human Reproduction*, 15:366-372
- Auger J, Eustache F** (2000) Standardisation de la classification morphologique des spermatozoïdes humains selon la méthode modifiée de David. *Andrologia*, 10:358-373
- Behre HM et al.** (2000) Diagnosis of male infertility and hypogonadism. In: Nieschlag E, Behre HM, eds. *Andrology, male reproductive health and dysfunction*. Springer: 92ff.
- Cooper TG et al.** (2010) World Health Organization reference values for human semen characteristics. *Human Reproduction Update*, 16:231-245
- Cross NL** (1995) Methods for evaluating the acrosomal status of human smear. In: Fenichel P, Parinaud J, eds. *Human smear acrosome reaction*. Paris, John Libbey Eurotext (Colloques INSERM): 277-285
- Kruger TF et al.** (1987) A quick, reliable staining technique for human smear morphology. *Archives of Andrology*, 18:275-277
- WHO Press** (2010) *Laboratory manual for the examination and processing of human semen*. 5th edition

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Referenzen

Andersen AG et al. (2000) High frequency of sub-optimal semen quality in an unselected population of young men. *Human Reproduction*, 15:366-372

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Behre HM et al. (2000) Diagnosis of male infertility and hypogonadism. In: Nieschlag E, Behre HM, eds. *Andrology, male reproductive health and dysfunction*. Springer: 92ff.

Cooper TG et al. (2010) World Health Organization reference values for human semen characteristics. *Human Reproduction Update*, 16:231-245

Cross NL (1995) Methods for evaluating the acrosomal status of human sperm. In: Fenichel P, Parinaud J, eds. *Human sperm acrosome reaction*. Paris, John Libbey Eurotext (Colloques INSERM): 277-285

Kruger TF et al. (1987) A quick, reliable staining technique for human sperm morphology. *Archives of Andrology*, 18:275-277

WHO Press, (2010) Laboratory manual for the examination and processing of human semen. 5th edition

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Behre HM et al. (2000) Diagnosis of male infertility and hypogonadism. In: Nieschlag E, Behre HM, eds. *Andrology, male reproductive health and dysfunction*. Springer. 92ff.

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Cross NL (1995) Methods for evaluating the acrosomal status of human smear. In: Fenichel P, Parinaud J, eds. *Human smear acrosome reaction*. Paris, John Libbey Eurotext (Colloques INSERM): 277-285

Kruger TF et al. (1987) A quick, reliable staining technique for human smear morphology. *Archives of Andrology*, 18:275-277

WHO Press (2010) *Laboratory manual for the examination and processing of human semen*. 5th edition