



Data sheet

(please read carefully)

ANTIGENES GmbH
Labordiagnostika

Römerstraße 20-22
58332 Schwelm
Tel.: +49 2336-9154958
Fax: +49 2336-9154957
Email: info@antigenes.de
Homepage: www.antigenes.de

Mites extract

Order-No.: AEA 5500	Lot: see Label	Quantity: 1 mg/ml
Source	Deutschland	
Raw material	Dust mites mix	
Consistance	Liquid (Tris-HCl buffer)/Cristalic (lyophilised)	
Additif preservation	0,1 % Sodium azide	
Date of manufacture	12.05.2021	
Expire	December 2025	
Proteine concentration analysis	determined by Bradford assay (BioRad) using BSA as standard	
Molecular weight determined	SDS-PAGE CE-Experion (BioRad)	
Size	1 mg/vial	
Storage	Shipped at 4°C, store at -18°C/8-15°C	

Application: for research and development uses, not as remedy, not for in vivo uses.

Foam or bubble development should be avoided as much as possible.

For optimal immunoassay use, we recommend performing a well-titration first to avoid false positive results.

The liquid antigen may be stored for a short time (a few days) at 8 °C. Unused stock solution must then be frozen at -18 ° C. Multiple freezing and thawing are not recommended because these may reduce significantly the antigen activity.

Schwelm, den 30.05.2021

This electronic document is valid without signature